



A Q U A S A F E®

## EVACUATION SYSTEM CLEANER

Bacterial-Based Cleaner is Safe for Equipment, the Environment, and Staff

- Bacterial-based cleaner provides residual overnight debris removal from cuspidor drains and suction systems
- Works under all conditions, even with cold water
- Neutral pH, safe for amalgam separators and traps
- AquaSafe® leaves a fresh, clean scent at traps and throughout the entire evacuation system
- Environmentally friendly - 6 non-pathogenic bacteria provide safe and natural cleaning power
- AquaSafe® ensures a level of protection that cannot be matched by other evacuation system cleaners that simply mask odor

AQSE - 64 oz. bottle	\$54.99 ea.
AQSE4 - Case of 4 bottles	\$199.95

\*Prices are subject to change.

### INCREDIBLY ECONOMICAL - 256 ONE-QUART TREATMENTS PER GALLON

#### Protection and Peace of Mind

An evacuation system cleaner is essential to safe and effective operation of every dental practice. Regular use of AquaSafe® reduces the presence of disease causing bacteria, offensive odors, and can reduce the chance of system failure resulting in a costly service calls and repair.

## DENTAL UNIT WATER TESTING PROGRAM

The EPA Recommends Quarterly Testing of Dental Unit Water in Each of Your Operatories

It is recommended that testing of dental unit water be conducted quarterly to ensure it meets CDC guidelines for bacteria levels. Centers for Disease Control and Prevention recommends that all dental unit water be below 500 colony forming units of heterotrophic water bacteria per ml of water. This easy-to-use kit provides in-office results and describes how to resolve elevated bacteria levels. In addition, the kit provides assistance with establishing protocol for continued compliance.

**Features include:**

- In-office reporting
- Helps to meet CDC Guidelines
- Protects patients
- Online Infection Control Course with 1 CEU\*\*

WTK-A (1) Unit	\$24.95
WTK-B (4) Units	\$89.95
WTK-C (8) Units	\$159.95
WTK-D (12) Units	\$199.95
WTK-E (24) Units	\$372.95

\*Prices are subject to change.

\*\* To redeem your online Infection Control Course, please email [hptc-orders@hptcinc.com](mailto:hptc-orders@hptcinc.com) and reference "Water Testing Training".

